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- Johnson, J. S., et al. (2019) Evaluation of 16S rRNA gene sequencing for species and strain-level microbiome analysis. *Nature Communications*. 10(1), 5029.
- Gehrig, J., et al. (2022) Finding the right fit: evaluation of short-read and long-read sequencing approaches to maximize the utility of clinical microbiome data. *Microbial Genomics*, 8(3), 10.1099/mgen.0.000794.
- PacBio datasets and NCBI PRJNA750084
- Overview – Sequel systems application options and sequencing recommendations. PacBio documentation.
- Procedure & checklist – Full-length 16S amplification, SMRTbell library preparation and sequencing. PacBio documentation.
- Procedure & checklist – Preparing whole genome and metagenome libraries using SMRTbell prep kit 3.0. PacBio documentation.
- Shoreline Biome Kits | Shoreline Biome SB Analyzer
- <https://github.com/xfengnfx/hifiasm-meta>
- Kolmogorov, M., et al. (2020) metaFlye: scalable long-read metagenome assembly using repeat graphs. *Nature methods*, 17(11), 1103–1110.
- Koren, S., et al. (2020) HiCanu: accurate assembly of segmental duplications, satellites, and allelic variants from high-fidelity long reads. *Genome Research*. 30, 1291-1305.
- Callahan, B. J., et al. (2019) High-throughput amplicon sequencing of the full-length 16S rRNA gene with single-nucleotide resolution. *Nucleic Acids Research*, 47(18), e103.
- <https://qiime2.org/>
- [https://github.com/LangilleLab/microbiome_helper/wiki/PacBio-CCS-Amplicon-SOP-v1-\(qiime2\)](https://github.com/LangilleLab/microbiome_helper/wiki/PacBio-CCS-Amplicon-SOP-v1-(qiime2))
- <https://www.onecodex.com/>
- <https://ezbiome.com/>
- <https://github.com/PacificBiosciences/pb-16S-analysis/wiki/Analyzing-PacBio-HiFi-Mock-Community-16S-Data-with-QIIME-2>
- Huson, D. H., et al. (2018) MEGAN-LR: new algorithms allow accurate binning and easy interactive exploration of metagenomic long reads and contigs. *Biology direct*, 13(1), 6.
- <https://github.com/PacificBiosciences/pb-metagenomics-tools>
- <https://bugseq.com/>

 Learn about metagenome sequencing:
pacb.com/microbial-pop

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Information in this document is subject to change without notice. PacBio assumes no responsibility for any errors or omissions in this document. Certain notices, terms, conditions and/or use restrictions may pertain to your use of PacBio products and/or third party products. Refer to the applicable PacBio terms and conditions of sale and to the applicable license terms at <http://www.pacb.com/legal-and-trademarks/terms-and-conditions-of-sale/>. PacBio, the PacBio logo, SMRT, SMRTbell, and Sequel are trademarks of PacBio. All other trademarks are the sole property of their respective owners.

© 2022 PacBio. All rights reserved. Research use only. Not for use in diagnostic procedures.
102-193-684 REV02 05MAY2022

PacBio